

Spectroscopic imaging system for high-throughput viability assessment of ovarian microtumors in a microfluidic system

Amélie St-Georges-Robillard^{1,2}, Mathieu Masse^{1,2}, Jennifer Kendall-Dupont^{2,3}, Mathias Strupler¹, Michael Jermyn^{4,1}, Anne-Marie Mes-Masson^{2,3,5}, Frédéric Leblond^{1,2}, Thomas Gervais^{1,2,3}

¹Department of Engineering Physics, Polytechnique Montréal, Canada;

²Centre hospitalier de l'Université de Montréal Research Center, Canada

³Institut du cancer de Montréal, Canada;

⁴Brain Tumour Research Centre, Montreal Neurological Institute and Hospital, Department of Neurology and Neurosurgery, McGill University, Canada

⁵Department of Medicine, Université de Montréal, Canada

ABSTRACT

There is a growing effort in the biomicrosystems community to develop a personalized treatment response assay for cancer patients using primary cells, patient-derived spheroids, or live tissues on-chip. Recently, our group has developed a technique to cut tumors in 350 μm -diameter microtissues and keep them alive on-chip, enabling multiplexed in vitro drug assays on primary tumor tissue. Two-photon microscopy, confocal microscopy and flow cytometry are the current standard to assay tissue chemosensitivity on-chip. While these techniques provide microscopic and molecular information, they are not adapted for high-throughput analysis of microtumors.

We present a spectroscopic imaging system that allows rapid quantitative measurements of multiple fluorescent viability markers simultaneously by using a liquid crystal tunable filter to record fluorescence and transmittance spectra. As a proof of concept, 16 microtissues composed of ovarian cancer cell line TOV112D were trapped in a microfluidic system, stained with a live cell marker (CellTrackerTM Green or Orange) and a dead cell marker (propidium iodide), and imaged. System response, background noise and autofluorescence were removed. Fluorescence intensity was normalized by the white-light spectrum at excitation to isolate the markers' intrinsic fluorescence from tissue absorption and scattering. Spectral un-mixing was applied to separate each fluorophore's contribution. We have demonstrated that rapid and simultaneous viability measurements on multiple microtissues can be achieved, which will have a significant impact to predict a tumor's response to multiple treatment options as well as in drug discovery to assess the potential of a drug candidate directly on human primary tissue.

This poster was presented during the:

SPIE Photonics West BIOS 2016

Sunday Poster Session

Sunday, February 14th, 2016

Paper Number: 9689-160; Time 5:30 PM - 7:00 PM